(12) INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(19) World Intellectual Property

Organization

International Bureau



(43) International Publication Date 12 October 2017 (12.10.2017)

- (51) International Patent Classification: *G01N33/50* (2006.01) *C12Q 1/37* (2006.01)
- (21) International Application Number: PCT/US20 17/026 172
 - 101/05201//0201/2
- (22) International Filing Date: 5 April 2017 (05.04.2017)
- (25) Filing Language: English
- (26) Publication Language: English
- (30) Priority Data: 62/3 18,659 5 April 2016 (05.04.2016) US
- (71) Applicant: UNIVERSITY OF PITTSBURGH OF THE COMMONWEALTH SYSTEM OF HIGHER EDUCATION [US/US]; 1st Floor Gardner Steel Conference Center, 130 Thackeray Avenue, Pittsburgh, PA 15260 (US).
- (72) Inventors: ACHARYA, Abhinav; 2428 Bear Run Dr., Pittsburgh, PA 15237 (US). LITTLE, Steven, R.; 9804 Braewick Drive, Allison Park, PA 15 101 (US). TARIN, Tatum, V.; 1350 Lucust St., Building C Suite G100A, Pittsburgh, PA 15219 (US).
- (74) Agent: RUPERT, Wayne; Klarquist Sparkman, LLP, One World Trade Center, Suite 1600, 12 1 SW Salmon Street, Portland, OR 97204 (US).

(10) International Publication Number WO 2017/176899 Al

- (81) Designated States (unless otherwise indicated, for every kind of national protection available): AE, AG, AL, AM, AO, AT, AU, AZ, BA, BB, BG, BH, BN, BR, BW, BY, BZ, CA, CH, CL, CN, CO, CR, CU, CZ, DE, DJ, DK, DM, DO, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN, HR, HU, ID, IL, IN, IR, IS, JP, KE, KG, KH, KN, KP, KR, KW, KZ, LA, LC, LK, LR, LS, LU, LY, MA, MD, ME, MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO, NZ, OM, PA, PE, PG, PH, PL, PT, QA, RO, RS, RU, RW, SA, SC, SD, SE, SG, SK, SL, SM, ST, SV, SY, TH, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, ZA, ZM, ZW.
- (84) Designated States (unless otherwise indicated, for every kind f regional protection available): ARIPO (BW, GH, GM, KE, LR, LS, MW, MZ, NA, RW, SD, SL, ST, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, RU, TJ, TM), European (AL, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV, MC, MK, MT, NL, NO, PL, PT, RO, RS, SE, SI, SK, SM, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, KM, ML, MR, NE, SN, TD, TG).

Published:

— with international search report (Art. 21(3))

(57) Abstract: A Fe(II)- and Ca⁺² - chelated alginate/gelatin conjugate.



ASSAY FOR DETECTION OF BLADDER OR PROSTATE CANCER

This application claims the benefit of U.S. Provisional Appl. No. 62/318,659, filed April 5, 2016, 5 which is incorporated herein by reference.

BACKGROUND

Currently bladder and prostate cancers are the most expensive malignancy because of the regular use of cystoscopy for the follow-up of these patients. In order to decrease the cost from the use of cystoscopy, several urine markers have been developed over the years for regular follow-up. Clinically available tests include NMP-22, urine cytology, FISH, BTA-STAT among others. All these tests are send-off lab tests and are not available for in-clinic use and therefore there is a need for development of technologies that can be used in-clinic in order to decrease the costs. Moreover, these techniques need to have similar sensitivity to existing tests to be able to reliably monitor bladder cancer progression.

- 15 Such techniques need to detect biomarkers of bladder and prostate cancers. Several biomarkers for these cancers have been studied. Among them up-regulation of MMP2 and MMP9 (72 kDa and 92 kDa type IV collagenases or gelatinase A and gelatinase B, respectively) has been shown to be involved in cancer metastasis and therefore, is an attractive target to study bladder cancer progression. Moreover, clinical detection of these enzymes in biological fluids can differentiate between aggressive and non-aggressive
- 20 disease states and be used to monitor disease progression. These MMPs are associated with malignant phenotype of cancer cells because of their ability to degrade type IV collagen, which is a major component of the basement membrane. Moreover, MMP2/9 have been shown to have a strong correlation to unfavorable prognostic factors of urothelial bladder cancer and tumor recurrence and therefore may be useful in the follow up of the patients.
- 25 Quantitative detection of MMP2/9 is currently performed using enzyme linked immunosorbent assays (ELISAs - detection limit of 10⁻¹² M), whereas detection of the active form of enzyme is typically performed using zymography, (detection limit of 30x10⁻¹⁵ moles). Although, ELISA is quantitative highly sensitive, it cannot determine the active form of the enzymes. On the other hand, zymography can detect the active form of the enzyme, but they are not very sensitive. Both the techniques require sophisticated
- 30 instruments and trained personnel, and therefore are difficult to implement in clinics. Moreover, it is important to differentiate between benign bladder/prostate disease and cancer, which has been strongly linked to the active form of MMPs in the urine.

SUMMARY

35 Disclosed herein is a Fe(II)- and Ca⁺²- chelated alginate/gelatin conjugate, which may take the form of particles.

- 1 -

Also disclosed herein is a composition that includes the conjugate. In particular, the composition is a bioassay that includes a plurality of conjugate particles. In certain embodiments, the particles may agglomerate.

5

Further disclosed herein is a kit comprising:

a composition comprising the conjugate; and

a composition comprising a Fenton reagent and at least one chromogenic substrate. In certain embodiments the conjugate composition and the Fenton reagent composition are provided in separate containers.

Also disclosed herein is method for detecting the presence of cancer in a subject, comprising: contacting a biological sample from the subject with the conjugate to form a sample/conjugate composition; and

15

10

contacting the sample/conjugate composition with a composition that comprises a Fenton reagent at least one chromogenic substrate resulting in an assay composition;

wherein detection of a color change in the assay composition indicates the presence of a matrix metalloproteinase in the biological sample. The presence of a matrix metalloproteinase in the biological sample may be an indicator of the presence of cancer.

20

25

Additionally disclosed herein is a method comprising:

contacting a biological sample from a subject with a gelatin-crosslinked aggregate of particles comprising a Fe(II)- and Ca⁺²- chelated alginate or with a Fe(II)-chelated alginate;

generating individual particles via gelatinase-mediated cleavage of the gelatin crosslinks between the particles, wherein the gelatinase is present in the biological sample; and

contacting the individual particles with an acid, H2O2, and at least one chromogenic substrate.

Further disclosed herein is a composition comprising a gelatin-crosslinked aggregate of particles comprising a Fe(II)- and Ca^{+2} - chelated alginate.

30

35

Also disclosed herein is a method comprising:

contacting an (i) alginic acid or alginate with (ii) a Fe(II) ion source and, optionally, (iii) a Ca^{+2} ion source resulting in particle comprising Fe(II)- chelated alginate or Fe(II)- and Ca^{+2} - chelated alginate;

contacting the particles with an activating agent resulting in activated particles; and

contacting the activated particles with a gelatinase-specific gelatin-type substrate resulting in a gelatin-crosslinked aggregate of particles.

- 2 -

The foregoing will become more apparent from the following detailed description, which proceeds with reference to the accompanying figures.

BRIEF DESCRIPTION OF THE DRAWINGS

- 5 FIG. 1 is a schematic showing the mechanism of detecting gelatinases via an assay for matrix metalloproteinases 2/9 (MMP2/9) (referred to herein as "Ammp"). Step 1. Fe(II), the catalyst for the Fenton' s reaction is "hidden" inside the polymer alginate, which chelates Fe(II), and the activity of the catalyst is temporarily turned OFF (Note Ca⁺² is also chelated to the alginate but is not shown in FIG. 1). Step 2. Nanoparticles are generated from this chelated Fe(II)-alginate complex with a large quantity of
- 10 catalyst Fe(II) per nanoparticle. Step 3. These nanoparticles are then precipitated by crosslinking to each other with gelatin as the linker. Step 4. When the precipitated clusters of nanoparticles are exposed to collagenase (active MMPs 2/9) during the assay (such as from urine samples from a bladder cancer patient), the clusters release individual nanoparticles that become re-suspended into the solution. The amount of released individual nanoparticles is proportional to the amount of active MMPs 2/9 Step 5. Next, the re-
- 15 suspended nanoparticle solution is separated and is added to the chromophore. Step 6. The previously "hidden" Fe (II) catalyst is activated by adding acid to remove the alginic acid chelate, and the assay very effectively generates visually observable color output proportional to the amount of active MMPs 2/9. Steps 4-6 are utilized to identify bladder cancer by detecting collagenase in patient's urine. FIGS. 2A-2C show the characterization of Ammps and size change in the presence of collagenase type IV.
- 20 (FIG. 2A) Alginate (Fe/Ca) particles that were generated via inverse emulsion of alginic acid and chelation of equimolar levels of ferrous chloride and calcium chloride were imaged using SEM. (FIG. 2B) Conjugation of alginate (Fe/Ca) particles with gelatin type A causes aggregation of particles as observed by SEM. (FIG. 2C). The average size of alginate (Fe/Ca) particles, alginate (Fe/Ca) gelatin conjugated particles and alginate (Fe/Ca) gelatin conjugated particles in the presence of collagenase was obtained
- using dynamic light scattering.
 FIGS. 3A-3C: Ammps have a limit of detection of 1 pg/mL in water and urine spiked with collagenase type IV. (FIG. 3A) Oxidation of TMB and IR783 dyes (i) in the presence or absence of alginate (Fe/Ca) particles (representative of n = 3). (ii) Color key for increasing oxidation of TMB and IR783 in the presence of Ammps with time is shown. (FIG. 3B) Water was spiked with collagenase type IV ranging from 1 µg/mL
- 30 to 0.1 pg/mL and Ammps was performed to detect the limit of detection in water (representative of n = 3). (FIG. 3C) Non-diseased urine was spiked with collagenase type IV ranging from 1 µg/mL to 0.1 pg/mL and Ammps was performed to detect the limit of detection in urine (representative of n = 3). FIG. 4: Ammps have a specificity toward bladder/prostate cancer. Ammps performed on bladder/prostate cancer patients were compared with patients with benign prostatic hyperplasia (BPH) to determine
- 35 specificity.

FIG. 5: Photographs of Ammps incubated in presence of a first chromogenic substrate, but no acid, and Ammps incubated in presence of a first chromogenic substrate and HC1.

- 3 -

FIG. 6: Photographs of Ammps incubated in presence of a first chromogenic substrate and a second chromogenic substrate.

FIGS. 7A-7C. Ammps have a limit of detection of 1 pg/mL in water and urine spiked with collagenase type IV. (FIG. 7A) Oxidation of IR783 dyes (i) in the presence or absence of alginate (Fe/Ca) particles

- 5 (representative of n = 3). (ii) Color key for increasing oxidation of IR783 in the presence of Ammps with time is shown. (FIG. 7B) Water was spiked with collagenase type IV ranging from $1 \mu g/mL$ to 0.1 pg/mL and Ammps was performed to detect the limit of detection in water (representative of n = 3). (FIG. 7C) Non-diseased urine was spiked with collagenase type IV ranging from $1 \mu g/mL$ to 0.1 pg/mL and Ammps was performed to detect the limit of detection in urine (representative of n = 3).
- FIG. 8. Ammps identifies patients with bladder cancer. Ammps was performed in a blinded fashion on patients' urine samples and divided into positive (orange-yellow color) and negative (green-purple color) for bladder cancer. Representative images of Ammps performed on patients' urine samples are shown. Code for patients is shown on top and the results are shown on the bottom. * After performing cystoscopy, it was observed that the false positive case had a benign tumor mass.
- 15 FIG. 9: Gelatin-crosslinked-alginate particles induce color change via Fenton's reaction. Gelatin-crosslinked-alginate particles are able to oxidize IR783 and change their color from green to yellow via Fenton's reaction (representative of n = 2) when 6N HC1 was added. Moreover, addition of 6N HC1 in the absence of Fenton's reaction changes the color of un-oxidized IR783 dye from green to purple. FIG. 10: Absorbance decrease of IR783 with time in the Fenton's reaction. The absorbance of IR783 at 560
- 20 nm decreases with time as the reactive oxygen species react generated in the Fenton's reaction degrade the IR783 molecule.

FIG. 11: Detection limit of collagenase type IV in urine and PBS. The activity of collagenase in 30% nondiseased urine sample (in PBS) and in PBS was determined using a substrate of collagenase DQ-FITC. The limit of detection for urine was determined to be 1 ng/mL and in PBS it was determined to be 2 ng/mL.

FIG. 12: The cut-off level demarcating negative and positive results from Ammps were determined by optimizing for the highest specificity and sensitivity values. Below the brown color values were chosen to be negative because of the highest sensitivity and specificity at this demarcation.

DETAILED DESCRIPTION

30 Assays

Bladder and prostate cancer are the most expensive malignancies to treat from diagnosis to death because of the use of regular and expensive send-off lab tests, which are challenging to implement in resource-poor settings. Disclosed herein is an assay for matrix metalloproteinases 2/9 (MMP 2/9) (assay is referred to herein as "Ammp" or "Ammps") that is capable of in-clinic monitoring of the progression of

35 bladder/prostate cancer in resource-poor settings. A visual signal (e.g., a change in color, for example, from purple to yellow) is generated by Ammp, which is based on the gelatinase activity of MMPs on the Ammp substrates. In certain embodiments, Ammp substrates are generated by crosslinking Fe (II) chelated alginate particles and gelatin, and the cleavage of gelatin by gelatinase frees the alginate Fe to generate visual signal via Fenton's reaction. In certain embodiments, Ammps have the limit of detection of 1 pg/mL of gelatinases activity in urine samples and have high specificity toward bladder/prostate cancer. Importantly, Ammps follow the ASSURED guidelines of the World Health Organization and therefore can have an impact in resource-poor settings (see Table 1 below).

5

Table	1
-------	---

Affordable	\$0.045 per test	
Sensitive	%	
Specific	%	
User-friendly	Utilizes urine (non-invasive).	
Rapid	Takes less than 30 min including analysis.	
Equipment free	Visual test. No equipment required.	
Delivered	No cold-chain required.	

The Ammp assay requires no instruments to measure the signal generated by the active form of MMPs, and therefore, can be utilized in clinics. The Ammp assay detects MMP2/9 in human urine samples within 30

10 min. and therefore, are useful for in-clinic monitoring or home monitoring of bladder cancer patients or prostate cancer patients, and can potentially be utilized to detect other types of cancers such as breast cancer, brain cancer, renal cancer, and liver cancer among others. Ammp utilize a reaction substrate of gelatinases such as MMP2/9. The Ammp reaction substrate includes gelatin type A crosslinked alginate containing Fe/Ca ions. In certain embodiments, the gelatin type A crosslinked alginate containing Fe/Ca ions is in the

15 form of particles. As used herein, "particle" encompasses spheres, capsules, particles, or rods. The amount of MMP2/9 activity may be compared to collagenase type IV obtained from bovine sources.

Ammps is an assay that is highly sensitive to collagenase activity (observed in urine of bladder cancer patients) because it is able to (a) temporarily turn OFF (hide) the catalytic activity of Fenton's reaction catalyst and (b) turn ON (make available) large quantities of the catalyst in the presence of active

- 20 collagenase. The mechanism by which Ammp detects MMP2/9 is shown in FIG. 1. Detection of bladder cancer is based on the collagenase activity of active MMPs 2/9 on gelatin-crosslinked-alginate nanoparticles. The Fe(II) catalyst is "hidden" in these alginate nanoparticles through chelation with alginates (steps 1-3). When patient urine is added that contains active MMPs 2/9, alginate nanoparticles are released and resuspended from the gelatin-crosslinked-alginate aggregates (step 4). The amount of re-suspended
- 25 nanoparticles generated is proportional to the amount of active MMPs 2/9. This solution of re-suspended particles is then added to the chromophore dye, and Fe (II) is released from the re-suspended alginate nanoparticles by addition of acid (steps 5-6). The chromophore dye becomes oxidized in the presence of Fe (II) via Fenton's reaction and generates a visually observable color (range of color: purple (un-oxidized dye) to yellow (completely oxidized dye), proportional to the amount of active MMPs 2/9).

In certain embodiments, alginate particles were generated using Ca and Fe as the crosslinking agents. These particles were then conjugated with gelatin type A. At room temperature, gelatin coagulates, causing the particles to aggregate. In the presence of gelatinases, gelatin is cleaved thereby generating free floating particles. These free floating particles are isolated, which then initiate Fenton's reaction to generate a color-differentiated product that is visually observed.

5

The reaction substrate was generated by first forming alginate particles from alginic acid using an inverse emulsion method, with Fe (II) (e.g., from ferrous (II) chloride) and Ca (e.g., from calcium chloride) as the metal ions that are sequestered (e.g., via chelation) by alginic acid. The functional groups of alginic acid that are involved in the chelation are -COOH and -OH groups. Alginates were used as a base material

- 10 to generate ultra-sensitive substrates for collagenase because of its unique property to chelate divalent cations such as Fe(II) (using -COOH and -OH groups) at a very high density. This is advantageous because a large amount of Fe (II) will lead to faster reaction and rapid generation of a large amount of reaction product (color). In certain embodiments, the cation density can range from 0.001 mg to 10 mg, whereas the particle density can range from 2 particle to 1,000,000,000 particles.
- 15 The alginate particles may be formed by contacting an alginic acid with a Fe(II) ion source and, optionally, a Ca+2 ion source. Illustrative Fe(II) ion sources include iron oxide, iron chloride, iron phosphate, iron (II) sulfate, ammonium iron sulfate, dichlorotetrakis(pyridine)iron, iron bromide, iron citrate, iron fluoride, iron iodide, iron molybdate, iron nitrate, iron oxalate, iron perchlorate, iron phosphate, iron pyrophosphate, iron tetrafluoroborate, potassium hexacyanoferrate(II), among others in different
- 20 hydration states. Illustrative Ca⁺² ion sources include calcium chloride, calcium carbonate, calcium acetate, calcium citrate, calcium sulfate, calcium pyrophosphate, calcium propionate, calcium phosphate, calcium oxide, calcium oxalate, calcium lactate, calcium hydroxide, calcium gluconate, calcium glubionate, among others in different hydration states. Alginate particles can be made with only chelated Fe(II), but the yield may be too low. Ca^{+2} provides better binding to alginate than Fe, and thus is used in conjunction with Fe(II). 25 However, the Ca⁺² does not participate in the Fenton's reaction.

All forms of alginic acid or alginate can be used to generate the particles. These include alginates with medium or low viscosity and alginic acid at different forms of purification from algae. Moreover, alginic acid can be in the form of salts such as sodium among others that can also be used to generate alginate particles. Alginic acid provides a dual advantage of the capability (1) of generating particles and (2) chelating Fe and Ca ions.

30

In certain embodiments, the inverse emulsion may be a water-in-oil emulsion. An aqueous solution of alginic acid is prepared. This solution of alginic acid (which constitutes the water phase) is added to an oil phase. Illustrative oils for the oil phase include alkanes or substituted alkanes such as octane, isooctane, dichloromethane, ethyl acetate, dimethylformamide, hexanes, ether, or chloroform, aromatics such as

toluene, and combinations thereof. The Fe(II) ion source and the Ca^{+2} ion source are added to the resulting 35 alginic acid/oil mixture thus forming a water-in-oil emulsion. For example, the composition of the dispersed phase may be water (100 uL to 100 mL), a water soluble detergent that can be Triton x-100, Tween 20,

- 6 -

5

PCT/US2017/026172

Tween 80, or Pluronics among others (100 uL to 100 mL of 0.1- 10% wt/wt), alginic acid (O.Img to 100 mg), and the Ca and Fe ions (0.1 mg to 100 mg). The continuous phase may include the oil (such as octane (ImL to 1000 mL)) and oil soluble detergents such as Brij30, Brij 35, Brij 56, Brij 58, or Span 80 among others (100 uL to 100 mL of 0.1- 10% wt/wt). Next isopropanol (1ml to 10000 mL) is added to cure the particles. Other curing agents such as various alkanols such as methanol, ethanol among others can be used as well. Next the particles thus formed are precipitated by centrifuging in a centrifuge for certain amount of

time (1,2,3,4,5,10,20,30 minutes) at (1,5,10,20,30,40,50,100 xG force) to isolate the particles. These particles are then freeze dried and used.

In certain embodiments, the resulting alginate (Fe/Ca) particles have an average particle size of at least 800 nm, more particularly at least 1000 nm. In certain embodiments, the resulting particles have an average particle size of not greater than 5000 nm, more particularly not greater than 1000 nm.


The particles may be activated to provide functional groups on the particles for conjugating with a gelatinase-specific gelatin-type substrate. Particle activation may be accomplished by contacting the particles with an activating agent such as glutaraldehyde, paraformaldehyde, formaldehyde, Ethyl-3-(3-

- 15 dimethylaminopropyl)-carbodiimide), Sulfo-N-hydroxysulfosuccinimide, N-hydroxysulfosuccinimide, Azabenzotriazole- 1-yl)-N,N,N,N' -tetramethyluronium hexafluorophosphate, (2-(1H-benzotriazol- 1-yl)-1,1,3,3-tetramethyluronium hexafluorophosphate), and 0-(Benzotriazol- 1-yl)-N,N,N',N'-tetramethyluronium tetrafluoroborate among others.
- The activated particles may be conjugated to gelatinase-specific gelatin-type substrate by reacting or 20 contacting the particles with a gelatinase-specific gelatin-type substrate. Gelatinase-specific gelatin-type substrates include, but not limited to, materials that can be converted into gelation (e.g., collagen) and material that constitute a part of the extracellular matrix (e.g, fibronectin). Illustrative gelatinase-specific gelatin-type substrates include gelatin A, gelatin B, mixed gelatin, collagen type I, III, IV, V, VII and X, fibronectin, laminin, aggrecan, link protein, elastin, vitronectin, tenascinp, SPARC, decorinr, myelin basic
- 25 protein, alPls, al-antichymotrypins, IL-lbeta, proTNF- a, IGFBP-3, substance P, casein, and similar materials. In certain embodiments, gelatin type A (0.001 mg to 10 mg) was used to conjugate to alginate particles (0.01 mg to 100 mg). The functional group of alginate that is conjugated is -COOH and -OH groups, whereas the functional group of -NH2 and -COOH are conjugated on the gelatin. In certain embodiments, in order to generate nanoparticles containing Fe (II) that are sensitive to collagenase activity,
- 30 alginate particles may be crosslinked to each other with gelatin type A (a substrate of collagenase) as a linker. This crosslinking may be achieved by first activating the alginate particles with EDC/sulfo-NHS. Next, the EDC/sulfo-NHS is removed, and gelatin is added to the particles (incubated for 16 hours) to generate inter-crosslinked structures that would precipitate in solution.
- In certain embodiments, the resulting alginate (Fe/Ca) gelatin conjugated particles have an average 35 particle size of at least 0.005 mm, more particularly at least 0.01 mm. In certain embodiments, the resulting particles have an average particle size of not greater than 0.10 mm, more particularly not greater than 0.01 mm.

- 7 -

The alginate (Fe/Ca) - gelatin conjugated particles are contacted (e.g., mixed) with a biological sample (e.g., urine, blood, plasma or serum) from a subject. If gelatinase activity (e.g., MMP2 and/or MMP 9) is present in the biological sample, the enzyme cleaves the gelatin from the particles, thus freeing the alginate particles containing Fe (II) and Ca(II) ions.

- 5 At least a portion of the biological sample that includes the alginate (Fe/Ca) gelatin conjugated particles is contacted (e.g., mixed) with a Fenton reaction substrate composition. The Fenton reaction substrate composition includes a Fenton reagent and at least one chromogenic substrate that changes color in the presence of the Fenton reagent. In certain embodiments, the Fenton reaction substrate composition includes H202 (e.g., 0.0001 M to 8 M), HCl (e.g., 0.05 N to 6 N), at least one first chromogenic substrate that
- 10 changes color from transparent to yellow, and at least one second chromogenic substrate that changes color from purple to yellow. In certain embodiments, HBr or H2S could be substituted for HCl. Fenton's reaction generates radical oxygen species such as hydroxyl radicals and singlet radical oxygens that then oxidize the first chromogenic substrate (initially transparent) causing the change in its structure, generating a pi-pi conjugation structure thereby generating a color (purple, which in the presence of Ammp substrates turns
- 15 yellow). The second chromogenic substrate also gets oxidized by the radical oxygen species that causes it to change its color from purple to yellow. In certain embodiments, Fenton reaction substrate composition includes a Fenton reagent and at least one chromogenic substrate that changes it color from purple to yellow. In certain embodiments, Fenton reagent and at least one chromogenic substrate that changes it color from purple to yellow. In certain embodiments, Fenton reagent and a single chromogenic substrate, wherein the single chromogenic substrate is the second chromogenic substrate (e.g.,
- 20 the second chromogenic substrate as described below, or a similar second chromogenic substrate) Illustrative first chromogenic substrates include those having the structure of:




wherein each R1 is independently selected from -C³/₄, =C³/₄, ≡CH, -OH, or -N³/₄; and each R2 is
independently selected from -N³/₄, -OH, or -COOH. An illustrative first chromogenic substrate is 5,5'-Tetramethylbenzidine (TMB)


- 8 -


Illustrative second chromogenic substrates include 2-[2-[2-Chloro-3-[2-[1,3-dihydro-3,3-dimethyll-(4-sulfobutyl)-2//-indol-2-ylidene]-ethylidene]-l-cyclohexen-l-yl]-ethenyl]-3,3-dimethyl-l-(4-

sulfobutyl)-3//-indolium hydroxide, inner salt sodium salt (IR783), and compounds having a structure of:









wherein n = 1 to 7; and each R1 is independently an alkyl, substituted alkyl, aryl, or substituted aryl; 10 and each R2 forms part of a ring, particularly an aromatic ring such as pyrrole, indole or pyridine.

30

35

PCT/US2017/026172

The alkyl may be a branched or unbranched saturated hydrocarbon group of 1 to 24 carbon atoms, such as methyl, ethyl, re-propyl, isopropyl, re-butyl, isobutyl, f-butyl, pentyl, hexyl, heptyl, octyl, decyl, tetradecyl, hexadecyl, eicosyl, tetracosyl and the like. A "lower alkyl" group is a saturated branched or unbranched hydrocarbon having from 1 to 6 carbon atoms. Preferred alkyl groups have 1 to 4 carbon atoms.

- 5 Alkyl groups may be "substituted alkyls" wherein one or more hydrogen atoms are substituted with a substituent such as halogen, cycloalkyl, alkoxy, amino, hydroxyl, aryl, alkenyl, or carboxyl. For example, a lower alkyl or (Ci-Ce) alkyl can be methyl, ethyl, propyl, isopropyl, butyl, iso-butyl, sec-butyl, pentyl, 3-pentyl, or hexyl; (C3 -Ce)cycloalkyl can be cyclopropyl, cyclobutyl, cyclopentyl, or cyclohexyl; (C3-Ce)cycloalkyl can be cyclopropylmethyl, cyclobutylmethyl, cyclopentylmethyl,
- 10 cyclohexylmethyl, 2-cyclopropylethyl, 2-cyclobutylethyl, 2-cyclopentylethyl, or 2-cyclohexylethyl; (Ci-Ce)alkoxy can be methoxy, ethoxy, propoxy, isopropoxy, butoxy, iso-butoxy, sec-butoxy, pentoxy, 3-pentoxy, orhexyloxy; (C2-C6) alkenyl can be vinyl, allyl, 1-propenyl, 2-propenyl, 1-butenyl, 2-butenyl, 3-butenyl, 1,-pentenyl, 2-pentenyl, 3-pentenyl, 4-pentenyl, 1-hexenyl, 2-hexenyl, 3-hexenyl, 4-hexenyl, or 5-hexenyl; (C2-Ce)alkynyl can be ethynyl, 1-propynyl, 2-propynyl, 1-butynyl, 2-butynyl, 3-butynyl, 1-
- 15 pentynyl, 2-pentynyl, 3-pentynyl, 4-pentynyl, 1- hexynyl, 2-hexynyl, 3-hexynyl, 4-hexynyl, or 5-hexynyl; (Ci-Ce) alkanoyl can be acetyl, propanoyl or butanoyl; halo(Ci-C6)alkyl can be iodomethyl, bromomethyl, chloromethyl, fluoromethyl, trifiuoromethyl, 2-chloroethyl, 2-fluoroethyl, 2,2,2-trifluoroethyl, or pentafluoroethyl; hydroxy(Ci-C6)alkyl can be hydroxymethyl, 1-hydroxyethyl, 2-hydroxy ethyl, 1hydroxypropyl, 2-hydroxypropyl, 3-hydroxypropyl, 1-hydroxybutyl, 4-hydroxybutyl, 1-hydroxypentyl, 5-
- 20 hydroxypentyl, 1-hydroxyhexyl, or 6-hydroxyhexyl; (Ci-C6)alkoxycarbonyl can be methoxycarbonyl, ethoxycarbonyl, propoxycarbonyl, isopropoxycarbonyl, butoxycarbonyl, pentoxycarbonyl, or hexyloxycarbonyl; (Ci-Ce)alkylthio can be methylthio, ethylthio, propylthio, isopropylthio, butylthio, isobutylthio, pentylthio, or hexylthio; (C2-Ce)alkanoyloxy can be acetoxy, propanoyloxy, butanoyloxy, isobutanoyloxy, pentanoyloxy, or hexanoyloxy.
- 25 The aryl may be a monovalent unsaturated aromatic carbocyclic group having a single ring (e.g., phenyl) or multiple condensed rings (e.g., naphthyl or anthryl), which can optionally be unsubstituted or substituted.

Other chromogens can be used for the second chromogenic substrate such as ICG, IR780, cyanine 2, cyanine 2.5, cyanine 5, and cyanine 7. Other chromogens can be used for the first chromogenic substrate such as, paranitrophenol, (2,2'-Azinobis [3-ethylbenzothiazoline-6-sulfonic acid]-diammonium salt), (o-phenylenediamine dihydrochloride), (o-nitrophenyl - β -D-galactopyranoside), CN (4-chloro-1-naphthol) and DAB (3,3'-diaminobenzidine tetrahydrochloride) among others.

In certain embodiments, the color change that occurs over time in the presence of gelatinase is from green to blue to purple to orange to yellow. The color change can be visually detected by an un-aided human eye. The color change can also be read and recorded on a spectrophotometer.

Detecting the presence of, or an alteration in the amount of, one or more of the proteins using the methods disclosed herein indicates the prognosis or diagnosis of the subject, or indicates if a therapy is effective for treating a subject.

The term "subject" includes both human and non-human subjects, including birds and non-human 5 mammals, such as non-human primates, companion animals (such as dogs and cats), livestock (such as pigs, sheep, cows), as well as non-domesticated animals, such as the big cats. The term subject applies regardless of the stage in the organism's life-cycle. Thus, the term subject applies to an organism in utero or in ovo, depending on the organism (that is, whether the organism is a mammal or a bird, such as a domesticated or wild fowl).

10

Kits are also provided. The kit can include the alginate (Fe/Ca) - gelatin conjugated particles, and optionally the Fenton reaction substrate composition. The kit can include a container and a label or package insert on or associated with the container. Suitable containers include, for example, bottles, vials, syringes, etc. The containers may be formed from a variety of materials such as glass or plastic. The container typically holds a composition including one or more of the particles. In several embodiments the container

15 may have a sterile access port.

> A label or package insert indicates that the composition is of use for evaluating the presence or progression of cancer, particularly bladder cancer or prostate cancer in a subject. The label or package insert typically will further include instructions for use, such as particular assay conditions. The package insert typically includes instructions customarily included in commercial packages of products that contain

20 information about the indications, usage, contraindications and/or warnings concerning the use of such products. The instructional materials may be written, in an electronic form (such as a computer diskette or compact disk) or may be visual (such as video files). The kits may also include additional components to facilitate the particular application for which the kit is designed. Thus, for example, the kit may additionally contain means of detecting a color change (e.g. a color comparison guide). The kits may additionally

25 include buffers and other reagents routinely used for the practice of a particular method. Several embodiments are disclosed in the following numbered clauses:

1. A Fe(II)- and Ca⁺²- chelated alginate/gelatin conjugate.

30

2. The conjugate of clause 1, wherein the gelatin is gelatin A.

3. The conjugate of clause 1 or 2, wherein the conjugate is in the form of a particle.

4. The conjugate of any one of clauses 1 to 3, wherein the gelatin is a gelatinase-specific gelatintype substrate selected from gelatin A, gelatin B, mixed gelatin, collagen type I, III, IV, V, VII and X, 35 fibronectin, laminin, aggrecan, link protein, elastin, vitronectin, tenascinp, SPARC, decorinr, myelin basic protein, alPls, al-antichymotrypins, IL-lbeta, proTNF- a, IGFBP-3, substance P, or casein.

5. A composition comprising the conjugate of any one of clauses 1 to 4.

6. The composition of clause 5, wherein the composition includes a plurality of conjugate particles.

5

7. The composition of clause 6, wherein the particles have an average particle size of at least 100 nm, more particularly at least 800 nm.

8. The composition of clause 6 or 7, wherein the particles have an average particle size of notgreater than 1000 nm, more particularly not greater than 800 nm.

9. The composition of any one of clauses 5 to 8, wherein the composition is in the form of a bioassay.

15 10. A kit comprising:

the composition of any one of clauses 5 to 8; and a composition comprising a Fenton reagent and at least one chromogenic substrate.

The kit of claim 10, wherein the composition comprising a Fenton reagent and at least one
 chromogenic substrate comprises H2O2, HCl, at least one first chromogenic substrate, and at least one
 second chromogenic substrate.

12. A method for detecting the presence of cancer in a subject, comprising:

contacting a biological sample from the subject with the conjugate of any one of clauses 1 to 3 to form a sample/conjugate composition; and

contacting the sample/conjugate composition with a composition that comprises a Fenton reagent and at least one chromogenic substrate resulting in an assay composition;

wherein detection of a color change in the assay composition indicates the presence of a matrix metalloproteinase in the biological sample.

30

25

13. The method of clause 11, wherein the composition comprising a Fenton reagent and at least one chromogenic substrate comprises H2O2, HCl, at least one first chromogenic substrate, and at least one second chromogenic substrate.

35

14. The method of clause 12 or 13, wherein the conjugate is in the form of a particle.

15. The method of any one of clauses 12 to 14, wherein the matrix metalloproteinase is matrix metalloproteinase 2, matrix metalloproteinase 9, or both matrix metalloproteinase and matrix metalloproteinase 9.

5 16. The method of any one of clauses 12 to 15, wherein the cancer is prostate cancer or bladder cancer.

17. The method of any one of clauses 12 to 15, wherein the cancer is prostate cancer.

10

18. The method of any one of clauses 12 to 15, wherein the cancer is bladder cancer.

19. The method of any one of clauses 12 to 18, wherein the cancer can be distinguished from benign prostate hyperplasia.

15 20. The method of any one of clauses 12 to 18, wherein the cancer can be distinguished from benign conditions causing hematuria.

Example 1

Ammp substrate: An inverse emulsion of water in octane was utilized to generate particles. A 10 mg/mL
solution of alginic acid (Fisher Scientific, Pittsburgh, PA) in DI ³/₄ 0 was generated and used as the water phase. Span 80 (1.25 mL; Fisher Scientific) was dissolved in 75 mL of octane (Fisher Scientific, Pittsburgh, PA), and used as the oil phase. This solution of alginic acid was then added to the oil phase under stirring at 7,000 rpm (Silverson L4RT-A, East Longmeadow, MA) and allowed to mix for 3 min. 35 mL of 700 mM CaCl₂ solution and 35 mL of 700 mM ferrous (II) chloride tetrahydrate (Fisher Scientific, Pittsburgh PA)

- 25 solution made in DI ¾ 0 was then added under stirring at 7,000 rpm and allowed to mix for 3 min. 2propanol (100 mL; Fisher Scientific, Pittsburgh PA) was then added to the mixture in order to cure the particles, and allowed to mix for 3 min. During the process of curing the particles are stabilized and the detergent that covers the particles is removed. The particles obtained were centrifuged (Eppendorf) at 2000xGs for 2 min and the supernatant was discarded. The particles were then re-suspended in 2-propanol,
- 30 incubated at room temperature for 5 min and then centrifuged again at 2000xGs for 2 min. The alginate

particles were then lyophilized and used for further experiments. The resulting particles were imaged using a scanning electron microscope (SEM) and the average

size of particles was observed to be 900 nm using dynamic light scattering (DLS) (FIG. 2A).

10 mg of alginate particles were re-suspended in 0.9 mL of DIH $_2$ O and 0.1 mL of 25%

35 glutaraldehyde (Sigma-Aldrich, St. Louis, MO) was added to the solution. This solution was incubated at room temperature for 5 min. Next, the particles were washed 3 times to remove excess glutaraldehyde by centrifuging at 2000X for 5 min and re-suspending them in fresh 1 mL DIH ₂O. The approximate number of

particles was calculated from the average diameter of particles 800 nm and alginic acid density of 1.6 g/cm³. Two equivalents of gelatin type A (Bloom 225, Fisher Scientific, Pittsburgh PA; re-suspended at 1 mg/mL in DIH2O at 37 °C) was added to the particles. The reaction was carried out for 16 h under stirring at 37 °C. The particles thus formed were centrifuged washed 3 times at 2000X for 5 min using DIH 2O at 37 °C. These particles were then lyophilized and used as is.

5

The amount of gelatin per mg of particles was determined using BCA analysis to be 192+42 µg/mg. The resultant crosslinked particles were then imaged using SEM and it was observed that the particles were now agglomerated as compared to the non-conjugated particles (FIG. 2B). This agglomeration was further confirmed by sizing the particles using DLS, where the size of the aggregates was observed to be 2 - 6 μ m

- 10 (FIG 2C). Next, it was determined if the size of particles is modified in the presence of collagenases activity. 1 mg/mL collagenase type IV in 500 µL PBS was incubated with 50 µL of Ammp substrates (3 mg/mL) for 25 min and the size of the particles was measured using DLS at 25 min post incubation. It was observed that the size of the particles decreased from 2-6 µm to 500 nm.
- 15 Fenton's reaction substrate: Fenton's reaction substrate consisted of 100 µ^{*}_ 6N HCl (Fisher Scientific, Pittsburgh PA), 100 µL TMB substrate (R&D Systems, Minneapolis, MN), 100 µL H₂O ₂ (R&D Systems, Minneapolis, MN) and 10 µL of 1 mg/mL IR783 in DIH 20 (Sigma Aldrich, St. Louis, MO).
- Urine collection: All the protocols utilized here for collecting urine from patients were approved by the IRB 20 committee of the Human Research Protection Office of University of Pittsburgh. The urine was collected at UPMC Mercy Hospital and utilized on the same day.

In order to determine if the Ammp substrates perform Fenton's reaction and generate visual color, Ammp substrate 50 uL (3 mg/mL) in DIH 20 was incubated in the presence of TMB substrate and the change 25 in color was noted. It was observed that even after 5 min of incubation no color was generated (FIG. 5). Using these data we hypothesized that the Fe chelated by carboxyl and hydroxyl groups in the alginate particles were not able to participate in the Fenton's reaction. Therefore, 50 µī of 1 N HCl was added to the reaction mixture in order to release Fe from alginate particles and the reaction was carried out. It was observed that within 5 min a yellow color was generated which is typically observed due to oxidation of

- 30 TMB (FIG. 5). Soluble 0.1 mg/ml ferrous chloride incubated with the TMB and H_2O_2 was utilized as the positive control, which showed immediate generation of yellow color. Moreover, since, the Ammps need to be able to detect MMP2/9 in the urine samples (transparent to yellow colored), it is of great benefit if the starting color of the Ammp substrate was different than yellow. Therefore, we added another substrate of the Fenton's reaction, IR783, which turns from green to transparent in the presence of reactive oxygen
- species generated in the reaction (FIG. 6). Next the range of color change that occurred in the presence of 35 Ammp substrate and TMB + IR783 was tested. Figure 3A demonstrates that in the presence of Ammp

substrate the color changes from green - blue - purple - yellow with time. This change is depicted in the form of a gradient for reference (FIG. 3B).

Next, it was determined if Ammp could detect collagenase in phosphate buffered saline (PBS). 1 mL of PBS was spiked with collagenase type IV to generate $1 \mu g/ml$, 1 ng/ml, 100 pg/ml, 10 pg/ml, 1 pg/ml

- 5 and 0.1 pg/ml concentrations in 1.5 mL Eppendorf tubes. A 50 μL of Ammp substrate was added to the PBS and incubated at room temperature for 25 min. Next, 300 μL of the solution was retrieved without disturbing the Eppendorf tube and added to another tube containing HC1. IR783, TMB, and H2O2 were then added to this Eppendorf tube. The change in color was observed for 3 min and a photograph was obtained using a camera phone. It was observed that Ammps were able to change the color of the solution from
- 10 purple to yellow (purple = negative; yellow = positive) so as to differentiate between 1 pg/mL and 0 pg/mL of collagenase type IV in PBS (FIG. 3B). Similarly, Ammps was tested to determine whether it could detect collagenase type IV in human urine samples. Non-diseased urine sample was spiked with collagenase type IV to generate 1 µg/ml, 1 ng/ml, 100 pg/ml, 10 pg/ml, 1 pg/ml, 0.1 pg/mL and 0.01 pg/ml concentrations in 1.5 mL Eppendorf tubes. A 0.3 mL of the urine, 0.7 mL of 10X PBS and 50 µî₂ of Ammp substrate was
- 15 reacted for 25 min and 0.3 mL of the solution was retrieved without disturbing the Eppendorf tube and added to another tube containing HC1, after which IR783, TMB, and H₂O₂ were added. It was observed that Ammps were able to differentiate between 1 pg/mL and the negative control of 0 pg/mL of collagenase type IV with naked eye (FIG. 3C).
- In order to determine if the Ammps could be used for detecting collagenase in the urine samples of 20 bladder/prostate cancer patients, urine samples from a clinic were obtained and Ammps was performed on the samples in a blinded fashion. Patients' samples were coded with patient's initials and the date the samples were donated. A total of 26 patients were tested for the presence of collagenase type IV in the urine of patients. The urine was collected from patients and utilized fresh within 1 h. 50 μr⁻_L of Ammp substrate (3 mg/mL in PBS), 0.3 mL of the urine and 0.7 mL of 10X PBS was added in an Eppendorf tube and
- 25 allowed to react for 25 min. Next 0.3 mL of this reaction mixture was pipette out into another Eppendorf tube containing IR783, TMB, HC1 and H2O2. The color of the solution at the end of 3 min was noted and plotted on a graph y-axis representing different colors of the solution and the patient code on the x-axis (FIG. 4A). It was observed that a total of 14 patients were positively identified as being positive for bladder cancer using Ammps and 12 patients negative for bladder cancer using Ammps. There were 11 true
- 30 positives, 19 true negatives, 1 false positive and 0 false negative for bladder and prostate cancer. Also the Ammps were able to distinguish between patients with hematuria with bladder cancer and hematuria without bladder cancer with 100% specificity and sensitivity.

Example 2

35 *Ammps substrate:* An inverse emulsion of water in octane was utilized to generate nanoparticles. A 10 mg/mL solution of alginic acid (Fisher Scientific, Pittsburgh, PA) in DI H2O was generated and used as the water phase. Span 80 (1.25 mL; Fisher Scientific) was dissolved in 75 mL of octane (Fisher Scientific,

- 15 -

Pittsburgh, PA), and used as the oil phase. This solution of alginic acid was then added to the oil phase under stirring at 7,000 rpm (Silverson L4RT-A, East Longmeadow, MA) and allowed to mix for 3 min. Thirty-five (35) mL of 700 **mM** CaC12 solution and 35 mL of 700 mM ferrous (II) chloride tetrahydrate (Fisher Scientific, Pittsburgh PA) solution made in DI H20 was then added under stirring at 7,000 rpm and

- 5 allowed to mix for 3 min. One hundred (100) mL of 2-propanol (Fisher Scientific, Pittsburgh PA) was then added to the mixture in order to cure the particles, and allowed to mix for 3 min. The particles obtained were centrifuged (Eppendorf) at 2000xGs for 2 min and the supernatant was discarded. The particles were then re-suspended in 2-propanol, incubated at room temperature for 5 min and then centrifuged again at 2000xGs for 2 min. The alginate particles were then lyophilized and used for further experiments.
- 10 Ten (10) mg of alginate particles were resuspended in 0.9 mL of DIH2O and 0.1 mL of 20 mg EDC (1-ethyl-3-(3-dimethylaminopropyl)carbodiimide hydrochloride and 25 mg sulfo-NHS (N-hydroxysulfosuccinimide) (Sigma-Aldrich, St. Louis, MO) was added to the solution. This solution was incubated at room temperature for 25 min. Next, the particles were washed 3 times to remove excess EDC/sulfo-NHS by centrifuging at 2000X for 5 min and resuspending them in fresh 1 mL DIH2O. The approximate number of particles was
- 15 calculated from the average diameter of particles 700 nm and alginate density of 1.6 g/cm³, and 2 mole excess of gelatin type A to number of particles (Bloom 225, Fisher Scientific, Pittsburgh PA; re-suspended at 1 mg/mL in DIH2O at 37 °C) was added to the particles. The reaction was carried out for 16 h under stirring at 37 °C. The particles thus formed were centrifuged and washed 3 times at 2000X for 5 min using DIH2O maintained at 37 °C.
- 20 Fenton's reaction substrate: Fenton's reaction substrate consisted of 100 μL 6N HC1 (Fisher Scientific, Pittsburgh PA), 100 μL H₂0₂ (R&D Systems, Minneapolis, MN) and 10 μL of 5 mg/mL IR783 in DIH2O (Sigma Aldrich, St. Louis, MO). First the 0.3 mL of sample that had been incubated with Ammps for 25 min, was added to the 100 μ² 6N HC1 in an eppendorf tube. Next, 10 μ² of IR783 and 100 μ² H2O2 were added in this sequence. The reaction was carried out for 3 min and a camera was utilized to image the color
 25 generated
- 25 generated.

Urine collection and clinical data: All the protocols for urine collection were approved by the IRB committee of the Human Research Protection Office of University of Pittsburgh. The experiments were performed in a blinded fashion.

The results of Example 2 are shown in FIGS. 7-12.

and sensitivity were determined using the formulae:

30

The patients that were positive for cancer via conventional screening methods of cystoscopy, Gleason score and urine analysis were identified. Patients that had benign prostate hyperplasia were also identified using these tests. In order to determine specificity and sensitivity of Ammps, the false positive (FP) and false negative (FN) of the Ammps for detecting bladder cancer in patients was determined. FP and FN were determined by comparing the data obtained from Ammps with these conventional tests. Specificity

35

Number oftrue nægazives

Specificity = $\frac{1}{\text{number of true negatives} - \frac{1}{2} - \frac{1}{$

Number of true positives

Sensitivity = number of true positives + number of false negatives

Negative predictive value and positive predictive value are determined using the formula:

negative predictive value

Number of true aegat & 53/tained from cyst ssopj

number of true negatives obtained from cystoscopy + number of fsize negatives i from Ammps

positive predictive value

5

Number of true positives obtained from cystoscopy number of true positives obtained from cystoscopy + number of false positives from Ammps

- Ammps provided 100% sensitivity and 95% specificity in detecting cancer when the data was 10 compared to cystoscopy. The sensitivity and specificity of Ammps is comparable to BTA-STAT, Y and Z. Therefore, Ammps can have an impact on generating an assay that can be performed in-clinic. In order for Ammps to be applied in resource -poor settings, they should be able to follow the ASSURED guidelines recommended by WHO. Therefore, we performed cost-analysis for performing Ammps on a single patient. It was determined that the cost to perform each assay was 4.5 cents, which is affordable in several resource
- 15 poor settings. On the other hand, cystoscopy costs \$350-\$3000 and BTA-STAT (a urine biomarker assay) costs \$158 to \$228. Moreover, the sensitivity and specificity determined for Ammps was comparable with several of the conventional screening tests. Ammps can also be combined with these screening tests to improve the sensitivity and specificity of the bladder cancer tests. Ammps are also user-friendly because they utilize urine samples, which are obtained via non-invasive procedures, and are comparable to BTA-
- 20 STAT among others. In addition, Ammps could identify bladder cancer rapidly in less than 30 min, which is another criteria recommended by WHO, and other biomarker tests take more than 30 min to perform. For example, cystoscopy requires 1-2 h, urine cytology requires 2-3 days, and BTA-STAT requires 50 min. Ammps are performed in an instrument-free manner, which also brings down the total cost required to perform the test. Lastly, the reagents of Ammps can be stored in dry form at room temperature, and thus
 25 follows the delivered ordering of ASSURED.
- 25 follow the delivered criteria of ASSURED.

In view of the many possible embodiments to which the principles of the disclosed invention may be applied, it should be recognized that the illustrated embodiments are only preferred examples and should not be taken as limiting the scope of the invention.

30

What is claimed is:

1. A Fe(II)- and Ca⁺²- chelated alginate/gelatin conjugate.

5

2. The conjugate of claim 1, wherein the gelatin is gelatin A.

3. The conjugate of claim 1 or 2, wherein the conjugate is in the form of a particle.

4. The conjugate of any one of claims 1 to 3, wherein the gelatin is a gelatinase-specific gelatin-type 10 substrate selected from gelatin A, gelatin B, mixed gelatin, collagen type I, III, IV, V, VII and X, fibronectin, laminin, aggrecan, link protein, elastin, vitronectin, tenascinp, SPARC, decorinr, myelin basic protein, alPls, al-antichymotrypins, IL-lbeta, proTNF- a, IGFBP-3, substance P, or casein.

5. A composition comprising the conjugate of any one of claims 1 to 4.

15

6. The composition of claim 5, wherein the composition includes a plurality of conjugate particles.

7. The composition of claim 6, wherein the particles have an average particle size of at least 100 nm, more particularly at least 800 nm.

20

8. The composition of claim 6 or 7, wherein the particles have an average particle size of not greater than 1000 nm, more particularly not greater than 800 nm.

9. The composition of any one of claims 5 to 8, wherein the composition is in the form of a 25 bioassay.

10. A kit comprising: the composition of any one of claims 5 to 8; and a composition comprising a Fenton reagent and at least one chromogenic substrate.

30

11. The kit of claim 10, wherein the composition comprising a Fenton reagent and at least one chromogenic substrate comprises H2O2, HCl, at least one first chromogenic substrate, and at least one second chromogenic substrate.

35

12. A method for detecting the presence of cancer in a subject, comprising: contacting a biological sample from the subject with the conjugate of any one of claims 1 to 3 to

form a sample/conjugate composition; and

contacting the sample/conjugate composition with a composition that comprises a Fenton reagent and at least one chromogenic substrate resulting in an assay composition;

wherein detection of a color change in the assay composition indicates the presence of a matrix metalloproteinase in the biological sample.

5

13. The method of claim 12, wherein the composition comprising a Fenton reagent and at least one chromogenic substrate comprises H2O2, HCl, at least one first chromogenic substrate, and at least one second chromogenic substrate.

10

14. The method of claim 12 or 13, wherein the conjugate is in the form of a particle.

15. The method of any one of claims 12 to 14, wherein the matrix metalloproteinase is matrix metalloproteinase 2, matrix metalloproteinase 9, or both matrix metalloproteinase and matrix metalloproteinase 9.

15

16. The method of any one of claims 12 to 15, wherein the cancer is prostate cancer or bladder cancer.

17. The method of any one of claims 12 to 15, wherein the cancer is prostate cancer.

20

18. The method of any one of claims 12 to 15, wherein the cancer is bladder cancer.

19. The method of any one of claims 12 to 18, wherein the cancer can be distinguished from benign prostate hyperplasia.

25

20. The method of any one of claims 12 to 18, wherein the cancer can be distinguished from benign conditions causing hematuria.

21. A Fe(II)-chelated alginate/gelatin conjugate.

30

22. The conjugate of claim 21, wherein the gelatin is gelatin A.

23. The conjugate of claim 21 or 22, wherein the conjugate is in the form of a particle.

35 24. A composition comprising the conjugate of claim 23, wherein the composition includes a plurality of conjugate particles.

25. A kit comprising:

the composition of claim 24; and

a composition comprising a Fenton reagent and at least one chromogenic substrate.

5

26. A method for detecting the presence of cancer in a subject, comprising: contacting a biological sample from the subject with the conjugate of any one of claims 21 to 23 to form a sample/conjugate composition; and

contacting the sample/conjugate composition with a composition that comprises a Fenton reagent and at least one chromogenic substrate resulting in an assay composition;

10

wherein detection of a color change in the assay composition indicates the presence of a matrix metalloproteinase in the biological sample.


27. The method of claim 12, wherein the composition comprising a Fenton reagent and at least one chromogenic substrate comprises an acid and at least one chromogenic substrate.


15

28. The kit of claim 10, wherein the at least one chromogenic substrate changes color from purple to yellow.

29. The kit of claim 10, wherein the at least one chromogenic substrate is selected from 2-[2-[2-

20 Chloro-3-[2-[1,3-dihydro-3,3-dimethyl-1-(4-sulfobutyl)-2//-indol-2-ylidene]-ethylidene]-1-cyclohexen-1yl]-ethenyl]-3,3-dimethyl-1-(4-sulfobutyl)-3//-indolium hydroxide, inner salt sodium salt, or a compound having a structure of:





wherein n = 1 b 7; and each R1 is independently an alkyl, substituted alkyl, aryl, or substituted aryl;

5 and each R2 forms part of a ring, or

ICG, or

IR780, or

cyanine 2, or

cyanine 2.5, or

10 cyanine 5, or

cyanine 7.

30. The method of claim 12, wherein the at least one chromogenic substrate changes color from purple to yellow.


15

31. The method of claim 12, wherein the at least one chromogenic substrate is selected from 2-[2-[2-Chloro-34241,3-dihydro-3,3-dimethyl4^4-sulfobutyl)-2//4ndol-2-ylidene]-ethylidene]4-cyclohexen-l-

;or


; or

yl]-ethenyl]-3,3-dimethyl-l-(4-sulfobutyl)-3//-indolium hydroxide, inner salt sodium salt, or a compound having a structure of:



 R^{2} N^{+} R^{+} R^{+} R^{+} R^{+} R^{+} R^{+} R^{+} R^{+} R^{+} R^{+}

5



wherein n = 1 to 7; and each R1 is independently an alkyl, substituted alkyl, aryl, or substituted aryl; and each R2 forms part of a ring, or

10 ICG, or

;or


IR780, or cyanine 2, or cyanine 2.5, or cyanine 5, or

5 cyanine 7.


32. The method of claim 12, wherein the composition comprising a Fenton reagent and at least one chromogenic substrate includes only a single chromogenic substrate.


10

33. The method of claim 32, wherein the single chromogenic substrate is selected from 2-[2-[2-Chloro-3-[2-[1,3-dihydro-3,3-dimethyl-l-(4-sulfobutyl)-2//-indol-2-ylidene]-ethylidene]-l-cyclohexen-lyl]-ethenyl]-3,3-dimethyl-l-(4-sulfobutyl)-3//-indolium hydroxide, inner salt sodium salt, or a compound having a structure of:



15





wherein n = 1 to 7; and each R1 is independently an alkyl, substituted alkyl, aryl, or substituted aryl; and each R2 forms part of a ring, or

ICG, or

5 IR780, or

cyanine 2, or cyanine 2.5, or cyanine 5, or

cyanine 7.

10

34. A method for detecting the presence of cancer in a subject, comprising:

contacting a biological sample from the subject with the conjugate of any one of claims 1, 2, 3, 4,

21, 22 or 23 to form a sample/conjugate composition; and

contacting the sample/conjugate composition with an acid, H2O2, and at least one chromogenic

15 substrate resulting in an assay composition;

wherein detection of a color change in the assay composition indicates the presence of a matrix metalloproteinase in the biological sample.

35. The method of claim 34, wherein the conjugate is in the form of a particle.

20

36. The method of claim 34 or 35, wherein the at least one chromogenic substrate changes color from purple to yellow.

37. The method of any one of claims 34 to 36, wherein the acid is selected from HBr, 34S, HCl, or 25 a mixture thereof. 38. A method comprising:

contacting a biological sample from a subject with a gelatin-crosslinked aggregate of particles comprising a Fe(II)- and Ca⁺²⁻ chelated alginate or with a Fe(II)-chelated alginate;

5

generating individual particles via gelatinase-mediated cleavage of the gelatin crosslinks between the particles, wherein the gelatinase is present in the biological sample; and

contacting the individual particles with an acid, H2O2, and at least one chromogenic substrate.

39. The method of claim 38, wherein the gelatinase is matrix metalloproteinase 2, matrix10 metalloproteinase 9, or both matrix metalloproteinase and matrix metalloproteinase 9.

40. A composition comprising a gelatin-crosslinked aggregate of particles comprising a Fe(II)- and Ca^{+2} - chelated alginate.

15 41. A method comprising:


contacting an (i) alginic acid or alginate with (ii) a Fe(II) ion source and, optionally, (iii) a Ca⁺² ion source resulting in particle comprising Fe(II)- chelated alginate or Fe(II)- and Ca⁺²- chelated alginate; contacting the particles with an activating agent resulting in activated particles; and contacting the activated particles with a gelatinase-specific gelatin-type substrate resulting in a

20 gelatin-crosslinked aggregate of particles.


25

30

35



1/12






SUBSTITUTE SHEET (RULE 26)

SUBSTITUTE SHEET (RULE 26)


3/12




PCT/US2017/026172

PCT/US17/26172 01-05-2017








5/12


SUBSTITUTE SHEET (RULE 26)

FIG. 7A






7/12




WO 2017/176899


PCT/US2017/026172



8/12

SUBSTITUTE SHEET (RULE 26)





10/12

PCT/US17/26172 01-05-2017




FIG. 11

SUBSTITUTE SHEET (RULE 26)




FIG. 12

INTERNATIONAL SEARCH REPORT

International application No.

PCT/US2017/026172

A. CLASSIFICATION OF SUBJECT MATTER IPC (2017.01) GOIN 33/50, C12Q 1/37				
According to International Patent Classification (IPC) or to both national classification and IPC				
B. FIELD	DS SEARCHED			
Minimum documentation searched (classification system followed by classification symbols) IPC (2017.01) GoIN 33/50, C12Q 1/37				
Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched				
Electronic data base consulted during the international search (name of data base and, where practicable, search terms used) Databases consulted: Google Patents, CAPLUS, Google Scholar, PatBase Search terms used: Fe (II) and Ca+2 alginate/gelatin systems ,iron, calcium, conjugate composition, method for detecting bladder and prostate cancer				
C. DOCUMENTS CONSIDERED TO BE RELEVANT				
Category *	Citation of document, with indication, where ap	ppropriate, of the relevant passages	Relevant to claim No.	
E	YANG, Chen Y., et al. Fabrication of porous gelatin spinning process. Artificial Cells, Blood Substitutes, a 0 1 Jan 2009 (2009/01/01) abstract		1-9,21-24,40	
	XING, Qi, et al. Increasing mechanical strength of gelatin hydrogels by divalent metal ion removal. Scientific reports, 2014, 4:4706 01 Jun 2014 (2014/06/01) whole document		1,3-9,21-24,40	
A GUTIERREZ BANOS, J. L., et al NMP 22, BTA stat test and cytology in the diagnosis of bladder cancer: a comparative study. Urologia internationalis, 2001, 66.4: 185-190 0 1 Jan 2001 (2001/01/01) whole document		12-20,26,27,30-39, 4 1		
Further documents are listed in the continuation of Box C. See patent family annex.				
 * Special categories of cited documents: "A" document defining the general state of the art which is not considered to be of particular relevance "E" earlier application or patent but published on or after the international filing date "L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other "T" later document published after the international filing date or date and not in conflict with the application but cited to under the principle or theory underlying the invention "X" document of particular relevance; the claimed invention canno considered novel or cannot be considered to involve an invense step when the document is taken alone 		ation but cited to understand nvention uimed invention cannot be ed to involve an inventive		
special reason (as specified) "O" document referring to an oral disclosure, use, exhibition or other means "P" document published prior to the international filing date but later		Y document of particular relevance; the cla considered to involve an inventive step combined with one or more other such being obvious to a person skilled in the "&" document member of the same patent fa	when the document is documents, such combination art	
than the priority date claimed ac document memory of the same patent raining Daie of the actual completion of the international search Date of mailing of the international search report				
19 Jun 2017				
	ame and mailing address of the ISA: Authorized officer srael Patent Office VOLKOV Karina			
	Technology Park, Bldg.5, Malcha, Jerusalem, 9695101, IsraelFacsimile No. 972-2-5651616Telephone No. 972-2-5651777			