

Miltenyi MACSQuant
Current as of:

MACSQuantify Software v2.3
10/2/2015

Laser	Wave-length	new FL#	old FL#	Filter	Applicable Dyes and Stains
Violet	405	V-1	1	450/50	PacBlue, SytoxBlue, AF450, DCViolet, v450, eFluo450, Ho, BV421, VioBlue, VioletFluor
Violet	405	V-2	8	525/50	PacOrange, KrOrange, v500, AF430, VioGreen, BV510, PacGreen
Blue	488	B-1	2	525/50	FITC, GFP, AF488, CFSE, SytoxGreen, ALDH
Blue	488	B-2	3	585/40	PE, Dil, PI, Resazurin, PKH26
Blue	488	B-3	4	690/75	PerCP-Cy5.5, PE-Cy5.5, PI
Blue	488	B-4	5	750LP	PE-Cy7
Red	635	R-1	6	690/75	APC, AF647, TO-PRO3
Red	635	R-2	7	750LP	APC-Cy7, APC-H7

Service: 866-426-1194
Serial#: 2267